



ASAN PHARMACEUTICAL CO., LTD

163, Yeongcheon-ro, Hwaseong-si, Gyeonggi-do, 18462 Korea

TEL: +82-31-376-5990 FAX: +82-31-376-5993

www.asanpharm.com

Statement Letter

Dec.08, 2021

We, Asan Pharmaceutical Co., Ltd., are a manufacturer for the rapid diagnosis test kits for the products related to COVID-19 virus. We, hereby, would like to inform you of the result reading time of Asan Easy Test COVID-19 Ag.

We are pleased to confirm you that the reading time is 15minutes in accordance with Asan Easy Test COVID-19 Ag test insert. However, in the vast majority of cases, the results can be read in 3-5 minutes after the sample is applied to the test device (When the test line is shown within 1 minute, it can be read as positive result). Meanwhile, in rare cases reaction time may reach 15minutes. Therefore, in order to avoid any misinterpretation of the test results, we recommend checking a visual reaction of the test line appearance in 15 minutes after the sample is applied as well.

Sincerely,

A handwritten signature in blue ink that reads 'Yeom, Jeong Gyu'.

Yeom, Jeong Gyu

President

Asan Pharmaceutical Co., Ltd.